

**IMPORTANT INFORMATION  
CONCERNING YOUR NEW  
LAB-TECH GARMENTS**



Distributed by All-Spec Industries

5228 US HWY 421 N.  
Wilmington, NC 28401

Phone: 800-537-0351

Fax: 800-379-9903

web: [www.all-spec.com](http://www.all-spec.com)

email: [sales@allspec.com](mailto:sales@allspec.com)

**IMPORTANT NOTICE**

If your count does not agree with the packing list or if you are not fully satisfied, notify All-Spec immediately. All-Spec Industries cannot be held responsible for any discrepancies that are not reported within 7 days after receipt of order. Any garments that have been laundered, dry cleaned, altered, marked or embroidered cannot be returned.

**PRESCRIBED HOME CLEANING PROCESS**

Cold, cool or warm water washes with LIQUID laundry detergent (i.e. LIQUID Tide, LIQUID Dash, LIQUID Cheer, LIQUID Oxydol, etc.) Do not use laundry detergents, which contain added fabric softeners or bleach. In the case of severe soil or grease, stronger liquid detergents can be used, if they do not contain chlorine bleach.

**NO CHLORINE BLEACH  
NO FABRIC SOFTENER  
RECOMMEND HANG DRY**

Hang dry or tumble dry at lowest temperature setting on home dryer. Garment will dry very quickly. To ensure lasting performance, minimize the garment's exposure to excessive dryer heat.